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virus. HHV-6 hides in glandular tissue and in neural cells. The infection is widely distributed in general populations (3). IgG antibodies can be detected in about 80 % of people as a marker of previous infection. IgM antibodies are present both in acute primary infection and in reactivations. IgM cannot be considered as a marker of a recent infection (19). HHV-6 type B infection occurs in the first year of life. HHV-6 type A may occur later in life. Using the avidity of IgG both for HHV-6 type A and B can infections be found even after the first years of life. It appears that the spectrum of clinical symptoms associated with HHV-6 infection is enlarging and serological tests for both types of HHV-6 are needed to identify the primary and chronic infections associated with these viruses.

Principle of the Test

This indirect immunofluorescence test is a rapid, simple, sensitive method for detection of antibody to Human herpesvirus-6 both to type A and B as well as for differentiation of acute primary and chronic persistence infection. This fluorescence antibody assay utilizes the indirect method of fluorescence antibody staining first describe for more than 50 years ago. The test has been carried out in steps: first the human serum or cerebrospinal fluid is brought into contact with fixed infected cells. Antibody, if present in the test sample, will complex with the antigen in the cellular substrate. If the sample contains no antibody for this particular antigen, no complex is formed and all the serum components are washed away in the rinse cycle. The second step involves adding fluorescein-labeled anti-human antibody. If the antibody to HHV-6 is present (a positive reaction), bright apple-green fluorescence can be seen with the aid of a fluorescence microscope. In this assay two types of antigen coated slides are utilized, one for HHV-6 type A and one for HHV-6 type B.

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